

Alpha DESIGN® Liner Ordering Instructions

Want to order a DESIGN Liner? Here's what you need to know.

1. SELECT ONE OF THE FOLLOWING SHAPE CAPTURE OPTIONS

- Use **OMEGA® Tracer®**.
- Use the [Alpha DESIGN App](#) with an iPad and a Structure Sensor.
- Submit a cast or scan with a completed [online order form](#) from our website.

2. CASTING (IF APPLICABLE)

Please follow the [BK Casting Example](#) instructions on the website.

Please Note:



- For bilateral amputees, please use a separate order form for each residual limb to achieve the best socket/liner accuracy.
- When casting transtibial patients, please cast at least 4" to 5" above the knee.
- Please write the patient ID on the cast.

WillowWood ATTN: DESIGN Dept.
Ship casts to: 15441 Scioto Darby Road
Mt. Sterling, OH 43143 USA

3. ORDERING

Make sure to have the following information ready to enter into the software or the online order form:

- Measurements shown at right
- Purchase Order number
- WillowWood account number
- All patient and liner information
- Ship-to address
- .stl or .aop file (if applicable)
- **Note: Shape files can be uploaded to the online order form, or emailed to designlinerteam@willowwood.com**



After submitting the completed online order form, you will receive a summary of your order via email.
Please print the summary and include it in the shipping box if you are sending a cast.

Failure to complete and submit an online order form could delay your order.

WillowWood®

WillowWood Global LLC
15441 Scioto Darby Road
Mt. Sterling, OH 43143
phone 740.869.3377 / 800.848.4930
www.willowwood.com



WillowWood Global B.V.
Keizersgracht 62/64
1015 CS Amsterdam
The Netherlands
Patent www.willowwood.com/education-and-resources/patents

